

















	<p>Conformité Européenne. This product is compliant to the medical device directive 93/42/EEC (2007/47/EC). If followed by a number (Notified Body) it is a type approved device compliant to medical device directive 93/42/EEC (2007/47/EC) or in vitro diagnostic medical device 98/79/EC</p>
	<p>Use-by date. Do not use after date mentioned</p>
	<p>Caution, a specific warning or precaution associated with the device</p>
	<p>Manufacturer</p>
	<p>Contains or presence of natural rubber latex</p>
	<p>Batch code, Lot number</p>
	<p>Date of manufacture</p>
	<p>Catalogue number</p>
	<p>Sterilised using ethylene oxide</p>
	<p>Sterilised using irradiation</p>

	Do not re-sterilize
	Do not reuse. Re-use of the product can cause infection and other blood related transfusion diseases
	Temperature limitation for storage
	Keep away from (direct) sunlight
	Keep dry, avoid moisty environment
	Do not use is package is damaged
	In vitro diagnostic medical device
	Read instructions for use
	Suitable for use in the food industry